

PRODUCT SPECIFICATION

Nitrile Powdered Examination Gloves (Textured)

SECTION I: PRODUCT DESCRIPTION

1.1	Type	Nitrile Examination Glove, Powdered , Non-sterile
1.2	Material	100% Synthetic Nitrile Latex
1.3	Color (Pigment)	<ul style="list-style-type: none"> o Off White o Blue
1.4	Design and Feature	Ambidextrous, textured , beaded cuff
1.5	Powder	Inside absorbable cornstarch USPXX1
1.6	Storage Condition	The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight.
1.7	Shelf-Life	The gloves shall have shelf life of 5 years from the date of manufacture with the above storage condition.
1.8	Packing Style	100 pcs gloves x 10 dispensers x 1 carton
1.9	Size Marking	The size of gloves shall be marked in the check box on every carton with black ink.

SECTION II: PERFORMANCE REQUIREMENTS

(Sampling Plan – ISO 2859 Single Normal)

#	Characteristics	Inspection Level	Acceptable Quality Level	Reference Standard
2.1	Dimensions	S2	4.0	ASTM D6319-00a ^{e3}
2.2	Physical Properties	S2	4.0	ASTM D6319-00a ^{e3}
2.3	Freedom from Holes (Air Pump Test)	GI	1.5	ASTM D6319-00a ^{e3}
2.4	Visual Defects:			
(i)	Major Visual	GI	2.5	In-house practice
(ii)	Minor Visual		4.0	
2.5	Packaging Defects:			
(i)	Regulatory	GI	**	In-house practice
(ii)	Visual	GI	4.0	
(iii)	Critical (incl. Gloves Counting)	S2	4.0	
2.6	Powder Amount	N=2	-	ASTM D6319-00a ^{e3} ASTM D6124-06
2.7	Mix Size / Mix Glove / Mix Hand	Not Allowed		

** Unacceptable at any level

SECTION III: PERFORMANCE SPECIFICATION

3.1 Dimensions

Description	Size	Standard
Length (mm)	All Sizes	Min 240
Palm Width (mm)	XS	76 +/- 3
	S	84 +/- 3
	M	94 +/- 3
	L	105 +/- 3
	XL	113 +/- 3
Thickness (mm) *single wall	All Sizes	Finger: 0.11 +/- 0.02 (Typical value: 0.12 – 0.13) Palm: 0.09 +/- 0.02 (Typical value: 0.10 – 0.12)

3.2 Physical Properties

Description	Standard	
	Before Aging	After Aging
Elongation at break (%)	Min 500 (Typical value: 500 – 600)	Min 400 (Typical value: 400 – 500)
Tensile Strength (MPa)	Min 14 (Typical value: 14 – 16)	Min 14 (Typical value: 14 – 16)

3.3 Freedom from holes

The sample size and allowable number of non-conforming gloves in the samples shall be determined in accordance to Sampling Plan ISO 2859-1 Single Normal using inspection and acceptable quality level as stated in Section II: Performance Requirements.

3.4 Visual Defects

The sample size and allowable number of non-conforming gloves in the samples for both major and minor defects shall be determined in accordance to Sampling Plan ISO 2859-1 Single Normal using Inspection and acceptable quality level as stated in Section II: Performance Requirements.

3.5 Packaging Defects

The Sample size and allowable number of non-conforming in the samples for regulatory, visual and critical packaging defects shall be determined in accordance to Sampling Plan ISO 2859-1 Single Normal using inspection and acceptable quality level as stated in Section II: Performance Requirements (Gloves Counting=100 pcs by weight per Dispenser).

3.6 Powder Amount

Maximum 10 mg per dm²

Prepared by : Eva Vinoni Mustafa / QA Manager
 Authorized by : Noor 'Akilah Saidin / QA Deputy General Manager
 Date : 31st May 2010